

Combi Des Gel 85 %



PRODUCT TYPE	Disinfectant
APPLICATION	COMBI DES GEL 85% is used where a thorough and effective handdisinfection is required.
PROPERTIES	COMBI DES GEL 85% is an effective disinfectant due to content of ethanol, in the product. COMBI DES GEL 85% contains glycerine which softens the skin, the product is therefore very gentle towards the user.
STORAGE	Store in tightly closed original container in a well-ventilated place. Do not store near heat sources or exposed to high temperatures. Storage: 0 - 30°C Durability: 36 months
SAFETY	Please refer to the label and safety data sheet for information on safe use, handling and transport of the product.

INSTRUCTIONS AND DOSAGE

2-5 ml. COMBI DES GEL 85% is applied to dry and visually clean hands. (1-2 strokes, using metering pump).

Distribute COMBI DES GEL 85% between the fingers, the hand back, the hand palm and around the wrist.

Keep hands moist with COMBI DES GEL 85% for 30 seconds, and rub to dryness.

In the event of an outbreak of, for example Norovirus or Clostridium difficile it is recommended that hands should be washed with water and soap and afterwards disinfected with COMBI DES GEL 85%.

PRODUCTDATA

Colour	Colourless.
Physical state	Gel
Odour	With odor.
Bulk density	~0,85
pH Delivery state	7-9

Use biocides safely. Always read the label and product information before use.